

REAGENT PORTFOLIO

Biochemistry reagent menu

Pack volumes in mL · × indicates number of bottles · R1 / R2: reagent components

Liver function				20 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H511	ALB	R:40×3	Bromocresol green method	CE / NMPA
H512	TP	R:40×3	Biuret method	CE / NMPA
H514	ALP	R1:40×2; R2:20×1	NPP substrate-AMP buffer method	CE / NMPA
H515	ALT	R1:40×2; R2:20×1	Alanine substrate method	CE / NMPA
H516	AST	R1:40×2; R2:20×1	Aspartic acid substrate method	CE / NMPA
H518	γ-GT	R1:40×2; R2:20×1	GCANA substrate method	CE / NMPA
H519	DB	R1:40×2; R2:20×1	Diazonium salt method	NMPA
H5110	DB	R1:40×2; R2:20×1	Vanadate oxidation method	CE / NMPA
H5111	TB	R1:40×2; R2:20×1	Diazonium salt method	NMPA
H5112	TB	R1:40×2; R2:20×1	Vanadate oxidation method	CE / NMPA
H5113	CHE	R1:40×2; R2:20×1	Butylthiocholine substrate method	NMPA
H5114	AFU	R:40×3	CNPF substrate method	NMPA
H5115	TBA	R1:45×2; R2:30×1	Enzymatic cycling method	NMPA
H5116	ADA	R1:45×2; R2:30×1	Peroxidase method	NMPA
H5117	5'-NT	R1:45×2; R2:30×1	Peroxidase method	NMPA
H5118	mAST	R1:35×1; R2:14×1	Enzyme inhibition method	NMPA
H5119	CG	R1:40×2; R2:20×1	Latex immunoturbidimetric method	NMPA
H5120	PA	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H5121	MAO	R:40×3	Enzymatic cycling method	NMPA
H5122	GR	R1:45×2; R2:18×1	Glutathione substrate method	NMPA

Other assays				2 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H5103	SOD	R1:45×2; R2:30×1	Pyrogallol substrate method	NMPA
H5104	SA	R1:45×2; R2:30×1	Enzymatic method	NMPA

REAGENT PORTFOLIO

Biochemistry reagent menu

Pack volumes in mL · × indicates number of bottles · R1 / R2: reagent components

Blood lipids				15 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H521	CHOL	R:40×3	CHOD-PAP method	CE / NMPA
H522	CHOL	R1:40×2; R2:20×1	CHOD-PAP method	NMPA
H523	TG	R:40×3	Enzymatic method	CE / NMPA
H5215	TG	R1:40×2; R2:20×1	Enzymatic method	NMPA
H524	HDL-C	R1:45×2; R2:30×1	Direct method	CE / NMPA
H525	LDL-C	R1:45×2; R2:30×1	Direct / surfactant removal method	CE / NMPA
H526	NEFA	R1:45×2; R2:30×1	ACS-ACOD method	NMPA
H527	Lp(a)	R1:40×2; R2:20×1	Latex immunoturbidimetric method	CE / NMPA
H528	ApoA1	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H529	ApoB	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H5210	ApoE	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H5211	HCV	R1:40×2; R2:20×1	Enzymatic cycling method	NMPA
H5212	sdLDL-C	R1:45×2; R2:30×1	Direct clearance method	NMPA
H5213	Lp-PLA2	R1:40×2; R2:20×1	Enzymatic method	NMPA
H5214	Lp-PLA2	R1:40×2; R2:20×1	Latex immunoturbidimetric method	NMPA

Glucose metabolism				8 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H531	GLU	R:40×3	Glucose oxidase method	CE / NMPA
H532	GLU	R1:40×2; R2:20×1	Hexokinase method	CE / NMPA
H533	GLU	R1:40×2; R2:20×1	Glucose oxidase method	NMPA
H534	FMN	R1:45×2; R2:30×1	Tetrazolium blue method	NMPA
H535	HbA1c	R1:45×2; R2:30×1	Latex agglutination method	CE / NMPA
H536	GA	GA R1:40×2; R2:20×1 ALB R:40×3	Enzymatic method	NMPA
H537	ADPN	R1:40×2; R2:20×1	Latex immunoturbidimetric method	NMPA
H538	D3H	R1:40×3; R2:20×1	β-Hydroxybutyrate dehydrogenase method	NMPA

Pancreatic function				2 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H571	α-AMY	R1:40×2; R2:20×1	EPS-G method	CE / NMPA
H572	LPS	R1:40×2; R2:20×1	Enzymatic colorimetric method	CE / NMPA

REAGENT PORTFOLIO

Biochemistry reagent menu

Pack volumes in mL · × indicates number of bottles · R1 / R2: reagent components

Renal function				10 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H541	Urea	R1:40×2; R2:20×1	Urease–glutamate dehydrogenase method	CE / NMPA
H542	UA	R1:40×2; R2:20×1	Uricase–peroxidase coupling method	CE / NMPA
H543	Crea	R1:45×2; R2:30×1	Enzymatic method	CE / NMPA
H544	CYS-C	R1:40×2; R2:20×1	Latex immunoturbidimetric method	CE / NMPA
H545	β₂-MG	R1:40×2; R2:20×1	Latex immunoturbidimetric method	NMPA
H546	α₁-MG	R1:40×2; R2:20×1	Latex immunoturbidimetric method	NMPA
H547	NGAL	R1:40×2; R2:20×1	Latex immunoturbidimetric method	NMPA
H548	NAG	R1:45×2; R2:30×1	MPT substrate method	NMPA
H549	MALB	R1:40×2; R2:20×1	Immunoturbidimetric method	CE / NMPA
H5410	RBP	R1:45×2; R2:30×1	Latex immunoturbidimetric method	NMPA

Myocardial function				11 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H551	CK	R1:40×2; R2:20×1	Phosphocreatine substrate method	CE / NMPA
H552	CK-MB	R1:40×2; R2:20×1	Immuno-inhibition method	CE / NMPA
H553	LDH	R1:40×2; R2:20×1	Lactic acid substrate method	NMPA
H554	α-HBDH	R1:40×2; R2:20×1	α-Ketobutyric acid substrate method	NMPA
H555	LDH1	R1:40×2; R2:20×1	Lactic acid substrate method	NMPA
H556	Mb	R1:40×2; R2:20×1	Latex immunoturbidimetric method	NMPA
H557	cTnl	R1:40×2; R2:20×1	Latex immunoturbidimetric method	CE / NMPA
H558	H-FABP	R1:40×1; R2:40×1	Latex immunoturbidimetric method	NMPA
H559	IMA	R1:45×2; R2:30×1	Albumin–cobalt binding method	NMPA
H5510	CK-MB	R1:45×2; R2:30×1	Latex immunoturbidimetric method	NMPA
H5511	ACE	R:40×3	FAPGG substrate method	NMPA

Electrolytes & blood gas				11 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H561	Ca	R:40×3	Arsenazo III method	NMPA
H562	P	R:40×3	Phosphomolybdate method	NMPA
H563	CO₂	R1:40×1; R2:40×1	Phenol red colorimetric method	NMPA
H564	CO₂	R:40×3	PEPC enzyme method	NMPA
H565	Cl	R:40×3	Mercury thiocyanate method	NMPA
H566	Mg	R:40×3	Arsenazo III method	NMPA
H567	Fe	R1:40×2; R2:20×1	Ferrozine method	NMPA
H568	Cu	R1:40×2; R2:20×1	5-Br-PADAP method	NMPA
H569	Zn	R1:40×2; R2:20×1	5-Br-PADAP method	NMPA
H5610	K	R1:45×2; R2:30×1	Pyruvate kinase method	NMPA
H5611	Na	R1:45×2; R2:30×1	Galactosidase method	NMPA

REAGENT PORTFOLIO

Biochemistry reagent menu

Pack volumes in mL · × indicates number of bottles · R1 / R2: reagent components

Specific proteins				21 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H581	IgA	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H582	IgG	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H583	IgM	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H584	C3	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H585	C4	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H586	Fg	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H587	AAG	R1:40×2; R2:20×1	Immunoturbidimetric method	NMPA
H588	HPT	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H589	TRF	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H5810	AAT	R1:45×2; R2:30×1	Immunoturbidimetric method	NMPA
H5811	IgG4	R1:45×2; R2:30×1	Latex immunoturbidimetric method	NMPA
H5812	C1q	R1:40×2; R2:20×1	Immunoturbidimetric method	NMPA
H5813	ASO	R1:40×2; R2:20×1	Latex immunoturbidimetric method	NMPA
H5814	RF	R1:40×2; R2:20×1	Latex immunoturbidimetric method	CE / NMPA
H5815	CRP	R1:40×2; R2:20×1	Latex immunoturbidimetric method	NMPA
H5825	f-CRP	R1:40×2; R2:20×1	Latex-enhanced immunoturbidimetry	CE
H5817	Anti-CCP	R1:45×2; R2:30×1	Latex immunoturbidimetric method	NMPA
H5818	PCT	R1:45×2; R2:30×1	Latex immunoturbidimetric method	CE / NMPA
H5819	SAA	R1:45×2; R2:30×1	Latex immunoturbidimetric method	CE / NMPA
H5820	MMP-3	R1:45×2; R2:30×1	Latex immunoturbidimetric method	NMPA
H5824	KL-6	R1:45×2; R2:30×1	Latex immunoturbidimetric method	NMPA

Calibration & control				3 items
Cat. No.	Test item	Pack (mL)	Method	Certification
H591	CMCAL	5×4	—	CE / NMPA
H592	CMCON1	5×4	—	CE / NMPA
H593	CMCON2	5×4	—	CE / NMPA

Registration documentation is product- and destination-specific. Confirm the current assay documents and supplied pack configuration before ordering.