

Clear answers for modern animal care.

VETERINARY DIAGNOSTIC SOLUTIONS

About Wisenz

Wisenz Biotech Limited is an IVD commercialization and supply platform connecting international partners with selected diagnostic instruments, reagents and manufacturing resources from China.

Drawing on 10+ years of international IVD sales experience, we match products to market needs and coordinate product information, technical discussions, training arrangements and after-sales communication with suppliers.

For distributors, veterinary clinics and laboratories



VETERINARY · HEMATOLOGY

YS-H25

Multifunctional Blood Analyzer



A compact 25-parameter POCT hematology platform using AI microscopic imaging for WBC, RBC and PLT analysis, with colorimetric HGB measurement.

Only 10 uL

Venous or capillary blood

8-minute result

Fast single-sample workflow

No liquid circuit

Low maintenance; no waste liquid

AI morphology

Cellular recognition and imaging

POCT - FOUR SIMPLE STEPS**1. Collect**

Obtain 10 uL blood

2. Add sample

Use single-test pack

3. Run

Automated analysis

4. Report

Digital result output

25 reportable parameters**WBC group (11)**

WBC, NEU#, LYM#, MON#, EOS#, BAS#, NEU%, LYM%, MON%, EOS%, BAS%

RBC group (8)

RBC, HGB, HCT, MCV, MCH, MCHC, RDW-SD, RDW-CV

PLT group (6)

PLT, MPV, PDW, PCT, P-LCC, P-LCR

VETERINARY · HEMATOLOGY

YS-H5

5-Part Veterinary Hematology Analyzer



A high-throughput CBC+DIFF analyzer combining multi-angle laser scattering, flow cytometry, electrical impedance and cyanide-free hemoglobin measurement.

60 samples/hour

Efficient routine laboratory testing

25 parameters

CBC plus 5-part differential

<=20 uL

Small whole-blood sample

Connected workflow

Two-way LIS and large storage

ANALYSIS

CBC+DIFF veterinary configuration

SAMPLE MODES

Whole blood and pre-diluted

DISPLAY

10.4-inch touchscreen

STORAGE

100,000+ results

CBC+DIFF parameter menu

WBC differential (11)

WBC, NEU#, LYM#, MON#, EOS#, BAS#, NEU%, LYM%, MON%, EOS%, BAS%

RBC parameters (8)

RBC, HGB, HCT, MCV, MCH, MCHC, RDW-SD, RDW-CV

Platelet parameters (6)

PLT, MPV, PDW, PCT, P-LCC, P-LCR

YS-M4

Compact Real-Time PCR System



An ultra-compact four-channel real-time PCR system supporting both extraction-free testing and YS-AC20-assisted nucleic-acid extraction workflows.

750 g

120 x 120 mm footprint

Four channels

FAM, HEX, ROX and CY5

10 C/sec

Fast maximum ramp rate

Single sample

Focused on-demand testing

DIRECT / EXTRACTION-FREE

1. Collect & mix

Prefilled direct buffer

2. Dispense

25 uL into PCR reagent

3. Quick spin

Centrifuge 10 seconds

4. Amplify

Run on YS-M4

WITH NUCLEIC-ACID EXTRACTION

1. Collect

Validated specimen

2. Extract

Process on YS-AC20

3. Quick spin

Prepare reaction tube

4. Amplify

Run on YS-M4

Veterinary molecular assay menu

Core single-target assays

CPV; CDV; CAV (generic, type 1 or 2); CPIV; CCoV; FCoV; Toxoplasma gondii; Bordetella bronchiseptica; Mycoplasma spp.; MHV; Giardia; Cryptosporidium; Tritrichomonas foetus; Rickettsia spp.; Bartonella henselae; feline hemotropic Mycoplasma; canine Mycoplasma.

Feline panels

Respiratory 5-panel; diarrhea 5-panel; FCoV + FPV duplex; diarrhea 4-panel.

Canine & zoonotic panels

Canine respiratory 6-panel; canine diarrhea 4-panel; zoonotic 6-panel.

Assay formats, pack sizes, performance claims and regulatory availability vary by market.

VETERINARY · BIOCHEMISTRY

YS-C100

Automatic Veterinary Biochemistry Analyzer



A compact, one-button POCT chemistry analyzer using single-use reagent discs and approximately 150 uL of serum or plasma.

14 minutes

Automated sample-to-result testing

Single-use disc

Closed, convenient reagent format

150 uL

Small serum or plasma volume

4.5 kg

Compact and connected

POCT - ONE-BUTTON WORKFLOW**1. Prepare**

150 uL serum/plasma

2. Load disc

Add sample to reagent disc

3. Start

Insert and press one button

4. Report

Review or print results

Representative reagent-disc profiles**Full package**

Up to 23 items

Comprehensive exam

18 items

Preoperative comprehensive

11 items

Basic check

10 items

Electrolytes

8 items

Contact Wisenz for the current analyte composition of each disc, local availability and regulatory status.

VETERINARY · BIOCHEMISTRY

YS-C200

Veterinary Automatic Biochemistry Analyzer



A micro-volume automatic platform for routine veterinary chemistry, inflammation, coagulation and electrolyte testing, with an integrated touchscreen and printer.

120-140 uL

Whole blood, serum or plasma

Approx. 10 min

Efficient single-sample result

Multiple profiles

Chemistry and coagulation profiles

Integrated system

Android, printer and USB

POCT - SAMPLE TO RESULT

1. Collect

Whole blood/serum/plasma

2. Add sample

Load the selected test pack

3. Start

Automated incubation and test

4. Report

On-screen and printed result

Veterinary biochemical kit menu

Comprehensive (24)

ALB, TP, TB, A/G*, TG, Crea, BUN, B/C*, PHOS, Ca, CK, GLOB*, GLU, LPS, GGT, TBA, LDH, ALT, tCO₂, ALP, UA, TC, AST, AMY

Liver & kidney (18)

ALB, TP, TB, ALP, TC, A/G*, GLOB*, TBA, GGT, AST, ALT, K+, Crea, BUN, PHOS, Ca, B/C*, UA

Preoperative (14)

ALB, GLU, Crea, TP, BUN, ALP, ALT, A/G*, B/C*, GLOB*, CK, LDH, AST, TB

Canine inflammation (12)

c-CRP, ALB, TP, GLOB, A/G*, ALT, AMY, LPS, CK, Crea, BUN, B/C*

Feline inflammation (14)

f-SAA, ALB, TP, GLOB*, A/G*, GGT, ALT, ALP, TBA, AMY, CK, Crea, BUN, B/C*

Additional profiles

Liver function (11); kidney function (7); preoperative (10); coagulation: PT, APTT, Fib, TT; electrolyte profile (13).

VETERINARY · BIOCHEMISTRY

YS-CS10

Veterinary Semi-Auto Biochemistry Analyzer



A compact open-system chemistry analyzer for laboratories that require flexible reagent selection, programmable methods and adaptable veterinary workflows.

Open system

Configure tests with selected reagents

Flexible methods

End-point, kinetic and fixed-time

7 wavelengths

Plus one optional channel

200,000+ results

WLAN, USB and LIS support

DISPLAY

7-inch IPS touchscreen

INCUBATION

10 positions; 25/30/37 C

PROGRAMMING

300+ test programs

PRINTER

Built-in 57 mm thermal

Configurable veterinary test menu

Routine chemistry categories

Liver function | Kidney function | Metabolic screening | Lipids | Enzymes | Minerals and electrolytes | Inflammation-related chemistry

YS-CS10 is an open-system analyzer and therefore does not have one fixed proprietary assay menu. Available analytes depend on the selected veterinary reagents, calibrators and validated application settings.

HY-611 VET

Animal Urine Chemical Analyzer



Automated urine-strip reading for feline, canine and other animal samples, with cold-light optics and integrating-sphere measurement.

200 samples/hour

Efficient routine strip analysis

11 parameters

Matched HY-11A VET strips

100,000 results

Onboard data storage

Connected reporting

Built-in printer and LIS output

SIMPLE URINALYSIS WORKFLOW

1. Prepare

Collect urine as instructed

2. Apply sample

Prepare strip per instructions

3. Load & start

Place strip; activate reading

4. Review

Read, print or transfer results

11-parameter urine chemistry menu

Cells and infection-related markers

WBC - Leukocytes
NIT - Nitrite
BLD - Blood

Chemistry and metabolism

PRO - Protein; GLU - Glucose
KET - Ketone; BIL - Bilirubin
URO - Urobilinogen

Additional measurements

pH; SG - Specific gravity
VC - Ascorbic acid
Compatible strip: HY-11A VET

245 x 255 x 110 mm | 2.25 kg | AC 100-240 V, 50/60 Hz. Follow the applicable strip instructions for preparation and timing.

YS-AC20

Automated Nucleic Acid Extractor



A compact magnetic-bead system for veterinary DNA/RNA preparation, using prefilled reagents and guided touchscreen operation.

1-4 samples

Flexible low-throughput extraction

8-13 minutes

Reagent and protocol dependent

Prefilled reagents

Simplified reagent handling

4.3-inch touchscreen

Guided benchtop operation

SIMPLE SAMPLE PREPARATION

1. Prepare

Add sample to prefilled cuvette

2. Load

Position cuvettes and magnetic rod

3. Start

Automated binding, wash and elution

4. Collect

Prepare DNA/RNA for PCR

Matched extraction chemistry

MTQ036V veterinary reagent; TM-36 format, 36 tests per kit. Use the Animal program and the applicable kit instructions.

Specimens and selected specifications

Whole blood, serum, plasma, saliva, oral/nasal and fecal swabs, tissues and other kit-supported specimens.

Sample: 50-500 uL | Elution: 60-100 uL
190 x 244 x 252 mm | 4.7 kg | UV lamp

VETERINARY / MICROSCOPY

YS600 / YS800 / YS900

Biological Microscope Series



Three observation configurations for veterinary laboratory microscopy. All models can be paired with the HD3900 camera for live viewing and documentation.

Three configurations

Darkfield, phase contrast or dual mode

40x objective

Configuration-specific optical setup

Precision focusing

0.005 mm fine-focus precision

HD3900 imaging

2 MP; 1080p at 60 fps; HDMI

Choose the observation configuration

YS600 - Darkfield

40x plan achromatic objective.
Darkfield condenser.

YS800 - Phase contrast

40x phase-contrast objective.
Phase-contrast condenser.

YS900 - Dual mode

40x phase-contrast objective.
Pull-out darkfield / phase-contrast switching slide.

SHARED PLATFORM

30-degree trinocular head | 5 W LED illumination | Stage: 210 x 140 mm | X/Y travel: 75 x 50 mm
Interpupillary range: 48-75 mm | C-mount camera connection | HD3900 available with all three models

VETERINARY · RAPID TESTS

Veterinary Rapid Tests

81 colloidal-gold rapid tests organised across six species and application groups.



Canine · 31

- CDV Ag / Ab
- CPV Ag / Ab
- CCV Ag / Ab
- CAV-II Ag
- CAV-I Ag
- CAV Ab
- Canine Flu Ag
- CPV / CCV Ag
- CDV / CAV-II Ag
- CDV / Flu Ag
- Rabies Ag / Ab
- Heartworm Ag
- Rotavirus Ag
- Relaxin Ag
- Pancreatitis Ag
- Progesterone
- Ehrlichia Ab / Ag
- Lyme Ab
- Anaplasma Ab
- Babesia Ab
- Leptospira Ab
- Leishmania Ab
- IgE
- CRP
- TSH
- Cortisol

Feline · 17

- FPV Ag / Ab
- FCV Ag / Ab
- FHV Ag / Ab
- FCoV Ag
- FIPV Ag / Ab
- FeLV Ag
- FIV Ab
- FHV-I Ag
- Leishmania Ab
- Ehrlichia Ab / Ag
- Pancreatitis Ag
- SAA

Avian · 4

- Avian Influenza Ag
- H5 Subtype Ag
- Newcastle Disease
- Newcastle Disease Virus Ag

Bovine · 10

- Early Pregnancy Ag
- Progesterone
- IBR Ab
- Tuberculosis Ab
- BVDV Ag / Ab
- Herpesvirus Type I Ag
- Brucella Ab
- Coronavirus Ab
- FMD NSP Ab

Porcine · 9

- ASFV Ag
- CSFV Ag / Ab
- PCV2
- FMDV Ab
- Pseudorabies Ab
- PRRSV Ab
- Parvovirus Ag
- Brucella Ab

General Veterinary · 10

- Cryptosporidium Ag
- Entamoeba histolytica Ag
- Brucella Ab
- Giardia lamblia Ag
- Toxoplasma gondii Ag / Ab
- D-Dimer
- Troponin I
- Serum Amyloid A-like protein
- Anthrax Ab

QUICK-TEST WORKFLOW

1. Collect

Use the specified specimen

2. Apply

Add sample and buffer

3. Wait

Develop for stated time

4. Read

Interpret control/test lines

CONSUMABLES

Ultrasonic lysis

Volumetric dispensing tubes and nucleic acid release reagent.

A consumable set for ultrasonic direct amplification. The sample stays in the original tube, lysis takes about a minute, and pre-aliquoted or lyophilized PCR reagent is dispensed into the reaction tube. The same release step is also described for animal and plant epidemic samples.

About 1 minute

Ultrasonic lysis on the original tube.

Within 30 minutes

One sample addition, then a result.

Over 99%

Agreement with magnetic-bead extraction.

Under 2 μ L

Dispense error, including repeats.

OPERATION

1 Sample

Add the swab sample to the volumetric dropper, or another recommended tube.

2 Lyse

Load the original tube. One-click start. About 60 seconds.

3 Dispense

Pre-aliquoted or lyophilized PCR reagent, 10-25 μ L, into the PCR tube.

SPECIFICATIONS

Dispense volumes: 5-10 μ L, 15 μ L, 20 μ L and 25 μ L.

PCR load volume: 10-25 μ L.

Release reagent: room-temperature storage, matched to the volumetric dropper.

Tube compatibility: about 90% of ultrasonic instruments, as stated by the supplier.

