

YD-H5

# 5-Part Hematology Analyzer

Fully automated CBC + DIFF for routine hematology

**60**

samples/hour

 **$\leq 20 \mu\text{L}$** 

sample volume

**25 + 4**

reported + study

**10.4"**

touch screen

A compact benchtop system combining three-angle laser scattering, flow cytometry and electrical impedance, with an integrated touch-screen interface.

## Three-angle laser scattering

Semiconductor laser technology supports leukocyte classification.

## Flexible sample modes

Venous whole blood, peripheral whole blood and prediluted peripheral blood.

## 3D cell analysis

Scatter plots and RBC / PLT histograms support result review.

## Connected workflow

100,000+ stored results and bidirectional LIS with HL7 support.

## YD-H5 / PRODUCT DETAILS

# Technical specifications

|                            |                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------|
| Analysis mode              | CBC + DIFF                                                                                                       |
| Throughput / sample volume | 60 samples/hour / $\leq 20 \mu\text{L}$ in CBC + DIFF mode                                                       |
| Measurement technology     | Semiconductor laser scattering; flow cytometry; 3D scatter analysis; electrical impedance; cyanide-free HGB      |
| Sample modes               | Venous whole blood; anticoagulated peripheral whole blood; prediluted peripheral blood                           |
| Display                    | 10.4-inch LCD touch screen; $800 \times 600$ resolution                                                          |
| Result graphics            | Two scatter plots; one 3D scatter plot (viewable as two 2D plots); WBC/BASO scatter plot; RBC and PLT histograms |
| Storage                    | 100,000+ records, including histograms, scatter plots and patient information                                    |
| Interfaces / LIS           | 1 network port; 2 USB ports; bidirectional LIS with HL7                                                          |
| Optional features          | External USB printer; Wi-Fi connection; built-in RFID reader                                                     |
| Operating environment      | 10-30 °C; 20-85% humidity; 70-106 kPa; altitude $\leq 2000$ m                                                    |
| Power supply               | AC 100-240 V; 50/60 Hz                                                                                           |
| Dimensions / weight        | $280 \times 470 \times 460$ mm (W $\times$ H $\times$ D) / approximately 22 kg                                   |

## Measurement parameters

25 reported parameters

WBC, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-SD, RDW-CV, PLT, MPV, PCT, PDW, P-LCR, P-LCC, NEU%, LYM%, MON%, EOS%, BAS%, NEU#, LYM#, MON#, EOS#, BAS#

4 study parameters ALY%, ALY#, LIC%, LIC#

## Performance

| Item | Linear range                  | Carryover    | Repeatability (CV) |
|------|-------------------------------|--------------|--------------------|
| WBC  | $0-100 \times 10^9/\text{L}$  | $\leq 0.5\%$ | $\leq 2.0\%$       |
| RBC  | $0-8 \times 10^{12}/\text{L}$ | $\leq 0.5\%$ | $\leq 1.5\%$       |
| HGB  | 0-250 g/L                     | $\leq 0.6\%$ | $\leq 1.5\%$       |
| PLT  | $0-1000 \times 10^9/\text{L}$ | $\leq 1.0\%$ | $\leq 4.0\%$       |

For professional use. Optional configurations are subject to order confirmation.